

1. 1, 2, 3, 4,

50PCS

|     |                                                                                                                       |
|-----|-----------------------------------------------------------------------------------------------------------------------|
|     |                                                                                                                       |
|     | S : 96 * 124cm M : 100 * 130cm L : 107 * 136cm XL : 110 * 140cm                                                       |
|     | PP                                                                                                                    |
|     | PP50g                                                                                                                 |
|     |                                                                                                                       |
| MOQ | 5,000PCS                                                                                                              |
|     |                                                                                                                       |
|     |                                                                                                                       |
|     | 1PC / polybag, 25pcs / CTN, 5PCS / polybag, 25pcs / CTN<br>1PCS / polybag, 50PCS / CTN, 10PCS / polybag, 100PCS / CTN |
|     | 1 * 40HQ 15~20                                                                                                        |
|     |                                                                                                                       |
|     | T / T, L / C, D / A, D / P,                                                                                           |
|     |                                                                                                                       |
|     | ISO13485. CE, FDA, Type5 & A 6                                                                                        |
|     |                                                                                                                       |

50PCS



orient Honest  
MEDICAL INTERNATIONAL

