

4 ply 4 ply 4 ply, 4 ply 4 ply,

4ply 4 ply 4 ply 4 ply 4 ply

4 ply 4 ply 4 ply 4 ply 4 ply

1) 4 ply : 4ply, PP 4 ply + 4 ply 4 ply + 4 ply 4 ply + PP 4 ply

2) 4 ply : 19.5X7.5cm

3) 4 ply : 4 ply, 4 ply 4 ply

4) 4 ply : 4 ply 4 ply 4 ply, 4 ply 4 ply, 4 ply 4 ply

5) 4 ply 4 ply : 4 ply 4 ply 4 ply

4 ply 4 ply 4 ply 4 ply 4 ply

50 4 ply / 4 ply, 4 ply 4 ply : 18 \* 10 \* 9cm

2000 4 ply / CTN, CTN 4 ply : 52 \* 38 \* 38cm

4 ply 4 ply 4 ply.

4ply 4 ply 4 ply 4 ply 4 ply

MOQ	10 4 ply
4 ply	4 ply 4 ply
4 ply 4 ply	T / T, L / C, D / A, D / P, 4 ply 4 ply, 4 ply 4 ply, 4 ply
4 ply 4 ply	1,000,000 4 ply / 4 ply, 1 * 40 HQ 4 ply 4 ply 15 4 ply 4 ply 4 ply
4 ply 4 ply	4 ply 4 ply 4 ply 4 ply 4 ply
4 ply	4 ply 4 ply, 2 4 ply 4 ply 4 ply 4 ply

4 ply 4 ply 4 ply 4 ply 4 ply

1) 3ply 4 ply 2 4 ply 4 ply 4 ply

2) 2 4 ply 4 ply 1ply 4 ply 4 ply 4 ply

3) 4ply 4 ply 4 ply 4 ply

4 ply 4 ply 4 ply 4 ply 4 ply



orient Honest  
MEDICAL INTERNATIONAL

