

Product details

PP 500ml, 17.5cm x 9.5cm x 17.5cm

Material	PP 500ml + PP 500ml + PP 500ml
Weight (g)	50g
Dimensions (cm)	17.5 x 9.5 x 17.5
Volume (ml)	* 9.5CM 17.5, 500ml
Usage	APRI / 17.5 / 9.5 / 17.5 / 500ml / 17.5 / 9.5 / 17.5
Quantity (PC)	50 PC / 17.5, 2000 PC / CTN, 500ml
Storage	500ml / 17.5 / 9.5 / 17.5
Notes	500ml, 17.5, 9.5, 17.5

500ml, 17.5cm x 9.5cm x 17.5cm

1) 10 °C to 30 °C storage temperature, 500ml / 17.5cm x 9.5cm x 17.5cm

2) 500ml, 17.5cm x 9.5cm x 17.5cm.

Product display



orient Honest
MEDICAL INTERNATIONAL



The regular shipment way



Our business process

