

**All isolation gowns are made of the highest quality materials, providing excellent protection and superior comfort:**

- Non-linting materials are smooth and comfortable to eliminate skin irritation
- Highly breathable construction layers remain cool, even during long cases
- Various Style can be customized according to your different requirement

<b>Fabric Material</b>	PP+PE material
<b>Style</b>	With ties on collar and waist, knitted cuffs
<b>Dimensions size</b>	X-3XL
<b>Fabric weight</b>	30-40g, could be as your request
<b>Color</b>	Green color
<b>Certificate</b>	ISO13485, FDA
<b>MOQ:</b>	10000 pcs
<b>Packing design</b>	can be as your design
<b>Packing</b>	1pc/plastic paper bag,50pcs/ctn, or as your requirement
<b>Capability</b>	60,000pcs per day
<b>Delivery time</b>	7-30 days
<b><u>Applications</u></b>	laboratory, clinic, food industry, dust-free workshop,hospital

**The picture of disposable Level 1 isolation gown with knitted cuffs medical manufacturer**







