

SMS 120 \* 140 cm, acn 35-45g, ISO13485, FDA, CE

- SMS 120 \* 140 cm, acn 35-45g, ISO13485, FDA, CE
- SMS 120 \* 140 cm, acn 35-45g, ISO13485, FDA, CE
- SMS 120 \* 140 cm, acn 35-45g, ISO13485, FDA, CE


Material	SMS
Size	120 * 140 cm, acn 35-45g, ISO13485, FDA, CE
Weight	35-45g, ISO13485, FDA, CE
Standards	ISO13485, FDA, CE
MOQ	10000 PC
Unit	1pc / plastic bag, 50pcs / ctn, 60,000pcs
Lead time	7-30 days
Payment	T/T, L/C, D/P, D/A, etc.

1 pc / plastic bag, 50 pcs / ctn, 60,000 pcs





**Velcro  
back neck**



**Round collar**



**Tie on waist**



**Knitted cuff**

