

· 1, [] [] : OH-3000

· 2.Material : PP, PP + PE, SMS, Microporous

· 3. 🖂 : 🖂 , 🖂 , 🖂 ...

· 4 \square \square : 50pcs / box, 100pcs / ctn. \square \square \square \square \square \square

· 5, □□□ : ISO13485, FDA

- -00000 000 00 000 0000 00 0000 00000.
- $\hbox{-0000 00 00 00 0000 0000 0000}.$
- -000 0000 00 00 000 00 000 00 0 0 0

	PP, PP + PE, SMS, □□ □ □□□□□
	X-5XL
	30-60g, [][] [][]
	ISO13485, FDA
MOQ:	5000 PC
	50 [] / [][], 2000 [] / ctn, [][] [][] [][]
	□ □ 60,000pcs
	7-30 [