000 0 000 00 00 :

100pcs/ bag, 2000pcs/ ctn. | | | | | : 36*29*54cm., 1*40Container 240,0000pcs | | | | | | | | | | | |

	PP, PP+PE, CPE
	38x15cm, 40x15cm, 41x17cm
	100pcs/bag, 2000pcs/ctn
00 000	
	50000 PC
	3 0 000 00 000 00 0 0 0 0 0 0 0 0 0 0 0
	T/T, PayPal, Western Union, L/C, D/A, D/P, Escrow,
	□□ : 600,000pcs/day, 1*40hq□ 30 □ □□□ □□ □ □□□□.
FOB 🔲	





