

PP 40 \* 50 cm 20-60g, 1000 PC / bag, 60,000 PC 7 ~ 30 days

- 100% virgin PP resin.
- FDA approved for food contact.
- High strength and low cost.

Material	PP (Polypropylene)
Color	White / Customizable
Size (s)	40 * 50 cm (16 * 20 inches)
Weight (w)	20-60g, 1000 PC / bag
Thickness	0.2mm - 0.3mm
Standards	ISO13485, FDA
MOQ	10000 PC
Lead Time	7 - 30 days
Packaging	100pcs / bag, 1000 PC / bag, 60,000 PC / bag
Applications	Medical, Industrial, Food Packaging

PP 40 \* 50 cm 20-60g, 1000 PC / bag, 60,000 PC 7 ~ 30 days





