

□□□ □ □□ □□ □ □ □□□ □ □ □□□ □ □□ □□

□□□ □ □□ □□□ □□ □□

- 1) □□□ □ □□ □□ □□ □□□ □ □, □ □□ □□□ □□□ □□ □□□□□ □ □ □□□ □ □□.
- 2) □□□ □□□□ 400kg □ □□ □□□ □□.
- 3) □ □□ 6□□ □□ □□ □□, □□□ □ □□□ □□□ □□□ □□□ □ □□□□ □□ □□□.

□ □ □ □□□ □ □□ □□ □□

□□	□
□	□□
□□□	□□□ □□ □□ □□ □□□ □□□ □ □□ □□ □□,
□□	195cm * 85cm (□□□□□□ □□□ □□ □□□□ □ □□)
□□	□□, □□□□ □□, □□□
moq	5000 pc
□□	□
□□ □□	t/t, L/C, d/A, d/P, □□□ □□□, □□□□, □□□
□□ □□	10, 000pcs/□, 1 * 40 HQ □□□□□ 15 □□□ □□ □ □□
□□ □□	□ □□□ □, □□ □□ □□□□□ 30days
□□ □□ □□	□: 1 □□ pc □□ □□□ □□ □□, 50 pc □□ Ctn □□ □□ □□: 1 □□ pc □□ □ □, 50 pc □□ ctn □□□ □□□ □□

□ □□□ □□□ □ □□ □ □ □□□ □ □□□ □



orient Honest
MEDICAL INTERNATIONAL

