

Good quality Disposable non woven medical warming blanket non woven moving blanket

Disposable non woven medical blanket is made of non woven and Polyetser in side.

Provides warmth and comfort during first aid emergencies

Widely used for single-patient use

Perfect for Emergency Use,first aid kits, urgent care centers, hospitals, and other medical care facilities

This Material Is Latex Free

Material:

100% polyester inside,outside PP nonwoven

Feature:Disposable, Waterproof

Type:flat sheet, Co-Friendly Blankets

Shape:rectangular

Pattern:Striped

Use:Hotel, Disposable, Hospital

is\_customized:Yes

Weight:220-700g, 220-700g

Season:FOUR SEASONS

Product name:disposable medical blanket

Keywords:disposable

Color:blue,white or customized

Fabric:100% Polyetser[]PP non woven

Usage:Airline,hospital,nursing home

Logo:Accept Customized Logo

Design:Accept Custom Designs

Package:1PC/poly Bag



110 X 190cm  
200g  
羽绒被