

[KF94 4](#)
[KF94](#)

[KF94 4](#)
[KF94](#)

	1 : pp () 2 : 3 : 4 : pp
	4Ply
s IZE	20*8 CM
FABRIC W 8	30GSM + 25GSM + 25GSM+25GSM, C BE
COLOR	
	CE, CE, ISO13485, FDA
MOQ :	20000pcs
	.
	5 0 PCS/B OX, 200 0 PC S/CTN,
	PC
	7-30
	, , , ,

[KF94 4](#)
[KF94](#)





