

00 0 000 00 3 000 000 00 00-00

00 0000 00000 000 000 000, 00 00 000

00 00 000 00 00 000000 00

00000 00, 00 00 00 00

000

00 30 00 000-00 000 00 00

000 00	PP + PE, PP 00 PE 00
000	000 000, 000 00
00 00	120 * 140cm, 130 * 150cm
00	00, 000, 000 00 000 0000
000	CE, ISO13485, FDA
MOQ :	100,000pcs
00 000	000 00000 0 0 0000
00	10pcs / bag, 1000pcs / ctn, 00 000 00
00	00 50,000pcs
00 00	5-15 0
00	00, 000, 000, 00 00, 00000 000

00 3 000 00 000 00 00 00



