

1) 100% 3-ply 18+20+20G

- 1) 100% 3-ply 18+20+20G
- 2) 50 pieces per box, 40 boxes per CTN

2) 100% 3-ply 25+25+25G

- 1) 100% 3-ply 25+25+25G
- 2) 50 pieces per box, 40 boxes per CTN
- 3) 100% 3-ply 25+25+25G
- 4) 50 pieces per box, 40 boxes per CTN
- 5) BFE 99 %, CE / ISO13485, ISO II / ISO IIR

ackaging

17.5 * 9.5CM, 3ply, 18 + 20 + 20G	50 / box, 40 boxes / CTN	51 * 19 * 41cm	1 * 40 HQ (PCS) 3,450,000 PC
17.5 * 9.5CM, 3ply, 25 + 25 + 25G	50 / box, 40 boxes / CTN	51 * 37 * 30cm	2,400,000 PC

100% 3-ply

100%	PP + 3-ply + PP
100%	100% 3-ply
100%	EA rloop
100%	* 9.5CM 17.5, 50 pieces per box
100%	50 / box, 2000 / CTN, 100% 3-ply
100%	50 pieces per box, 40 boxes per CTN
100%	50 pieces per box, 40 boxes per CTN

100% 3-ply?

- 1) 100% 100% 100% 100% 100% 100% 100% 100%, 100%.
- 2) 100% 100%, CE, ISO 100%.
- 3) 100% : 1000,000pcs / 1, 1 40HQ 10-15 100% 100% 100% *;
- 24 100% 4) OEM 100% 100%, 100%.

100% 100% 100% 100%





orient Honest
MEDICAL INTERNATIONAL

