

TYPE IIR 3 plys

1) 2) 3) 4)

- 1) 2) 3) 4)
- 2) PP
- 3) ()
- 4)


Table

Item	Description
Material	PP
Model	TYPE IIR, N95, N99
MOQ	50
Options	T / T, L / C, D / A, D / P, , ,
Quantity	1,000,000 / , 1 * 40 HQ 15
Weight	, ,
Packing	50 / 2000 / CTN, ,

Table

- 1) 2) 3)
- 2)
- 3)

Table





orient Honest
MEDICAL INTERNATIONAL

