

Product details

Product description
Material: PE coated. With PE coating, it is more durable and resistant to stains.

Quantity: 1 / 50 / CTN
Weight: 5 / 50 / CTN, net weight.

Material: PP + PE45G

Weight: 45 GSM

Material: PE, PP, laboratorial, etc.

Material: S / M / L / XL / XXL / XXXL,

Material: PP, PE, etc.

Material: 1PCS / polybag, 50PCS / CTN

Material: 80% PE, PP, etc.

Material: CE, ISO13485, FDA

Product display



orient Honest
MEDICAL INTERNATIONAL






orient Honest
MEDICAL INTERNATIONAL



Our business process



The regular shipment way



FAQ

Q : 1. 貴司の製品はどの国で生産されていますか、どこですか？

A : はい、製品は中国で生産されています。2003年から生産を開始し、2013年にアメリカ市場に進出しました。

Q : 2. What is the price of your product?

回答：製品の価格はUSD18から20の範囲です。

Q : 3. 貴司の製品はISO, EN, CE, FDA (FDA)の規格を満たしていますか？

回答：はい、製品はFDA, CE, ISO13485の規格を満たしています。REACHの規制も満たしています。

Q : 4. Can I get a sample of your product? How much does it cost?

A : はい、サンプルは3~7日以内に無料で提供いたします。