



(  
1  
2  
5  
)  
□  
□  
□  
2  
c  
m  
□  
□  
□  
±

□□□ □□□ □□□ □□ :

□□ □□□ □□□ □□□ □□

□□ PP □□

□□ □ / □□□□ □□ / □□□□□ /, □□□ □□□□□□ □□ □□□ □ □□□□.

□□ S-6XL

□□ □□□, □□, □□, □□□, □□, □□, □□□, □□ □

□□ 5PCS / □□, 10 □□ / CTN

□□ □□□ □□ □□ □□ □ □□ □□□ □□□□□ □ □□□

MOQ (5)000 □

□□ 3 □ □□□ □□□□ □□□ □□ □ □ □□□□

□□ □□ T / T, □□□, □□□ □□□, L / C, D / A, D / P, □□□□,

□□ □□ □□□□□□ : 10,000 / □, 1 \* 40HQ 10-15 □ □□□ □□ □ □□□

FOB □□ □□ / □□ / □□ / □□□

□□□ □□□ □□□ □□ :



orient Honest  
MEDICAL INTERNATIONAL





orient Honest  
MEDICAL SUPPLY CO., LTD.

