

SMS Specifications :

Our SMS is made of high quality SMS fabric. It is soft, comfortable, breathable, and easy to care for. It is suitable for all seasons and occasions.

SMS Specifications :

Material	High quality SMS fabric, soft and comfortable, suitable for all seasons
Pattern	SMS, PP, PP / PE, SMMS, etc.
Color	White, Black, Grey, Blue, Green, Red, etc.
Size	S M L XL 2XL 3XL 4XL 5XL
Weight	150g, 180g, 200g, 220g, 240g, 260g, 280g
Quantity	5PCS / lot, 10 PCS / ctnor; 1PCS / lot, 50PCS / CTN
Lead time	3-5 days for sample, 15-20 days for bulk order
MOQ	5000 PCS
Payment	T/T, L/C, D/A, D/P, etc.
Shipping	10,000 / lot, 1 * 40HQ 10-15 days to most ports
FOB	China / Vietnam / India / etc.

SMS Features :

- 1) Soft and comfortable.
- 2) Breathable and easy to care for.
- 3) Suitable for all seasons and occasions.

Notes :

- 1, Price list.
- 2, Lead time for sample;
- 3, L / C, T / T, D / P etc;
4. FOB, CNF, CIF, EXW etc etc etc;
- 5, Payment terms etc etc etc.
- 6, Quantity etc etc. 3 PCS etc etc etc.
- 7, OEM services etc etc etc. etc, etc, etc etc etc etc;

- 8, 00 00. 00 00 00 0000 25 ~ 30 40HQ 0000 0.
- 9, 00 00 0000 000000. 00, 00, 00 00, 000 00.
- 10, 00 000 00 00. 000 000 00 00 0000000 000000 0000!

000 SMS 000 00 00 :





orient Honest
MEDICAL INTERNATIONAL

