

oversleeves :

- 1. Sewing, White, 20x40cm, 20x60cm
- 2. Sewing
- 3. White, 20x40cm, 20x60cm
- 4. Size : 20x40cm, 20x60cm

oversleeves, 20x40cm, 20x60cm

oversleeves :

- 1. 1 ~ 40 HO
- 2. CE ISO FDA
- 3. 24

oversleeves :

- 1. 3-5
- 2. 4
- 3. 100PCS / 10bags / CTN
- 4. 40 * 20 * 32cm

oversleeves :

Weight	35g
Material	21.5
Color	White, 20x40cm, 20x60cm
Quantity	100PCS / 10bags / CTN
Dimensions	40 * 20 * 32cm

oversleeves :



