

MOQ : 5,000pcs

MOQ	5,000pcs
Material	OPP
Size	* 120cm 90, 102 * 128cm
Weight	11kg / CTN
Lead time	15-20 days
Payment	T/T
Sample	1PC / OPP, 250PCS / CTN, 51 * 30 * 23cm / CTN
Weight	11kg / CTN

MOQ : 5,000pcs

MOQ : 5,000pcs

MOQ : 5,000pcs

MOQ : 5,000pcs

1. PE OPP, 100% PE OPP, 60g OPP, 100% PE OPP, 60g OPP, 100% PE OPP, 60g OPP.

2. OPP OPP, OPP OPP & A1 OPP necessary OPP OPP! OPP OPP FOLED, OPP OPP OPP OPP OPP OPP. OPP OPP OPP OPP OPP OPP.

MOQ : 5,000pcs



orient Honest
MEDICAL INTERNATIONAL





orient Honest
MEDICAL INTERNATIONAL

