

PP 100% 100% :

- 1. PP 100%
- 2. 20-60g / m2
- 3. / /
- 4. /
- 5. , ,
- 6. : 10,000 / , 1 * 40HQ 10-15
- 7. OEM 24 ,
- 8. CE, FDA ISO 9001

PP 100% 100% :

- 1. , ,
- 2. , , . , ,
- 3. ,

PP 100% 100% :

PP 100% 100% , , ; ; ; ; .

PP 100% 100% :

PP 100% 100%
 : PP
 : 20-50gsm
 : ;
 ;
 / / / /
 / ;
 ;

: , , laboratorial, workshop
 : , , ,
 MOQ 5000PCS
 : 5 10PCS / ; 50 100PCS / CTN;
 : , 80 % , , ,
 : T / T, , L / C, D / A, D / P, ,

: OEM /
 : CE, ISO13485
 : capacity : 10,000 / , 1 * 40HQ 10-15

PP 100% 100% 100% :



orient Honest
MEDICAL INTERNATIONAL

