

[headwear cross kf94 mask adults multicolored](#)

All masks are made of the highest quality materials, providing excellent protection and superior comfort:

- Non-linting materials are smooth and comfortable to eliminate skin irritation
- Highly breathable construction layers remain cool, even during long cases
- Super soft earloop makes comfortable wearing for long time

[KF94](#)

Material	Layer 1 : pp (non-woven) Layer 2 : Non-woven Layer 3 : Non-woven Layer 4 : pp (non-woven)
Material	4 Ply
Size	20*8 CM
FABRIC W 8	30GSM + 25GSM + 25GSM+25GSM, C BE
COLOR	White
Standards	CE, CE, ISO13485, FDA
MOQ :	20000pcs
Weight	10g
Quantity	50 PCS/B OX, 2000 PCS/CTN, 2000 PCS/CTN
Material	Non-woven PC
Lead time	7-30 days
Payment	T/T, L/C, D/P, D/A, Western Union, etc.

[KF94](#)



orient Honest
MEDICAL INTERNATIONAL



orient Honest
MEDICAL INTERNATIONAL





 orient Honest
MEDICAL INTERNATIONAL

