

□□□□ □□□□ □□□ □□□ □□□ □□ □□ :

- 1. □□ : PP50g
- 2. □□ : □□,
- 3. □□□ : □□ □□ / □□, □□ □□
- 4. □□ : S-XXXL □□□□ □□
- 5. □□ : 5PCS / polybag□, 50PCS / □□□□ □□□□ □□□
- 6.1 * 20 "FCL : 50000PCS (□□ □□□□ □□ □□ □□□□ □□)
- 7. □□ □□ : 10, 000PCS
- 8. □□□□ : CE, FDA, ISO13485

□□□ □□□ □□□□ □□□ □□□ □□ :

- 1.Hand □□, □□, knittted □□□□ □ □□ □□.
- 2.Non-, □□□□, □□ □□ □□ □□□□□□, □□□□ □□□□ □ □□□ □□, □ □□□, □□□□□ □□□□.
- □□□□ □□ □□ 3.Against □.

□□ □□ □□ □□□□ □□ □□□□ □□ □□, □□ □□, □□ □□, □□ □□, □□ □□, □□ □ □□□□ □□ 4.Usually □□□□ □□□□ □□.

LAB COAT □□□□ □ :

□ □□ □□ □□ (□□□□ □□□□ □□□□ □□□ □□)								
□□□□	□□	□	□	□□	XXL	XXXL	XXXXL	XXXXXL
□□	(100)	(103)	(107)	(110)	(112)	(113)	(120)	(120)
□□ LENGTH	(132)	(135)	137	(140)	(142)	147	(150)	(158)

□□ □□□□ □□□□ □□□ □□□ □□ :

□□: □□□□ □□□□ □□□ □□□ □□, □□□□ □□ □□ □□□ □□, □□ □□□□ □□□ □□ □□ □□

□□: PP □ □

□□ : 50gsm

MOQ: 1000 pcs;
1000 pcs per lot;
1000 / 1000 / 1000 / 1000 pcs per lot
/ 10000;
1000 10000 pcs per lot

MOQ: 1000, 1000, laboratorial, workshop 1000

MOQ: 1000 pcs per lot per lot
5000PCS

MOQ: 5 100PCS / lot; 50 100PCS / CTN;

MOQ: 1000 pcs, 80 % 1000, 1000, 1000, 1000 10000 pcs

MOQ: T / T, 1000, 1000 1000, L / C, D / A, D / P, 10000,

MOQ: OEM 1000 pcs per lot / 1000

MOQ: CE, ISO13485

MOQ: capacity : 10,000 / lot, 1 * 40HQ 10-15 1000 pcs per lot

MOQ 1000 1000 1000 1000 1000 : 1000



orient Honest
MEDICAL INTERNATIONAL





orient Honest
MEDICAL INTERNATIONAL

