

1. □□ : PP50g 2. □□ : □□,

3. 000 : 00 00 / 00, 00 00 4. 00 : S-XXXL 000 00

7. [] [] : 10, 000PCS 8. [] : CE, FDA, ISO13485

000 000 000 000 000 00 :

1.Hand $\square \square$, $\square \square$, knittted $\square \square \square$ \square $\square \square$.

___ __ __ __ 3.Against __.

LAB COAT [][] [] :

0 00 00 00 (000 000 0000 000 00)								
					XXL	XXXL	XXXXL	XXXXXL
	(100)	(103)	(107)	(110)	(112)	(113)	(120)	(120)
□□ LENGTH	(132)	(135)	137	(140)	(142)	147	(150)	(158)

00 0000 0000 000 000 00 :

DD:	
ΠΠ:	PP ∏

□□: 50gsm

/ 0000;

□□: □□, □□, laboratorial, workship □□

MOQ 5000PCS

□□: 5 □□ 10PCS / □□; 50 100PCS / CTN;

OD OOO : OEM O OOO OO / OOO

□**□:** CE, ISO13485

000 000 000 000 000 00 :



