

□□ □□□□□□ □□□ □□ □□ :

□□□□□□ □□□ □□ □□ PP □ □ □□□ □□, □□ □□, □□□□ □□

□ LABORATORY, □□□□ □□ □□ □□, particl □□□□□□□□ □□ □□□ □□□ □□□

□□.

□□ □□□□□□ □□□ □□ □□ :

□□□:	□□□□□□ □□□ □□ □□
□□:	PP □ □
□□ :	40gSM, □□ □□□ □□□ □□
□□□:	□□ □□ □□ □□; □□□lastic □□ □□ □□; □□ □□ / □ / □□□□ / □□ □□ □□ □ / □□□□□; □□ □□ □□ □□□ □□

□□:	□□, □□, laboratorial □□, workshop □□
□:	□□, □□□, □□ □□, □, □□
□OQ	5000PCS
□□ :	1 5PCS / □□; (50) □□ 100PCS/ CTN;
□□:	□□□ □□, 80 % □□, □□, □□, □ □□□ □□□□ □□□

□□ □□: T / T, L / C, D / P, □□□ □

□□ □□□ :	OEM □ □□□ □□ □□ / □□□
□□:	CE, ISO13485

□□□ □□ □□ :

- □□ □ 1.5, 000 □□ □□ □□□□□ □□□□ □□□ □□□
- 2. □□□ □□ □□□ □□ □□□ □□ □□ □□□ □□.
- 3. □□ □□□□ □□□ □□□ □□ □□□ □ □□□□□.

4. 00 0 000 0 00

0 00 00 000 000 00 :





orient Honest
MEDICAL INTERNATIONAL

