

□ □□ / □□ □□ :

□□ : PP. □□□

□□ : 40gsm □□ □□□ □□

□□ : □□ □□ □□ □□□ □□

□□ : □□□ / □□□ / □□ □

□□ □□□□ : □□, □□□, □□, □□ □

□□□ : □□ □□□ □□□□ 80 % □□□ □□, □□, □ □□□ □□ □□

□□ : □□

□□ □ □□ : □

□□:

CE, ISO□□ □□□□□ □□ □□ □□.

2.PP □□□ □□ □□□ □□.

□□□□□ □□□□□ 3.Factory □□ □□.

4.Perienced □□, □□ □□, □□□□ □□ □□.

5. 10 □ □□ □□.

6.OEM □□ □□, □□ □□, □□, □□, □□ □ □□ □

7. □□ □□ □□ □□□□ □□□□□.

□□ □ □□ :

1.Normal □□□ : 1pcs / □□, 5pcs / mid bag, 600pcs / carton.

2.MOQ : 10000pcs.

3.Port : □□□ / □□, □□.

4. □□ □□□ : □.

5. □□ □□ □□ : □□ □□ 3 □, 25 □ □□ □□.

□□	□□□, SBPP, SMS, Spunlace, □□, □□ □□□
□□	□□ □□
□□	25-50gsm.
□□	□□□, □□□ □, □□□ □□, □□ □□□ □, □□□, □□□□, □□□□, □□□, □□□
□□□ □□ □□	1pcs / □□, 5pcs / mid bag, 600pcs / carton.
□□ □□	T / T, PayPal, Western Union, L / C, D / A, D / P, □□□□
□□	□□, □□, □□, □□□ □, □□ □□, □□, □□, □□ □□, □□□ □□ □
□□ □□ □□	□□ / □□□□ / □□ / □□□□, □□□□ □□ □□

□□ □□□ □□





orient Honest
MEDICAL INTERNATIONAL









