

Packaging of face mask with logo:

Specification	Normal packing	Carton size	1*40 HQ container can hold (Pcs)
17.5*9.5cm, 3-ply, 18+20+20G	50 pcs/box, 40 boxes/ctn	52*19*42cm	3,280,000 Pcs
17.5*9.5cm, 3-ply, 25+25+25G	50 pcs/box, 40 boxes/ctn	51*37*30cm	2,400,000 Pcs

Remarks: all can be special customized according to your actual requirements

Specifications of face mask with logo:

Origin	Xiantao, Hubei, China
Material	PP non woven+filter+pp non woven (Very low resistance to breathing,free from foreign smell,without any stimulating feeling)
Weight	15+20+25g/20+20+25g/20+25+25g (comfortable,without glass fibres,latex free.)
Style	Elastic ear-loop/Tie on (latex free)
Size	17.5x9.5cm(adult)/14.5x9.5cm(children)/12x7cm(baby)
Color	white,green,blue,pink,yellow,pink,orange,purple,red(all colors can be available)
Packing design	All inner box and carton printing can be as your design
Carton size	52x38x30cm / 52x19x42cm
Applications	hospital,laboratory,food industry,dust-free workshop,electronic manufactures and so on

Detail pictures of face mask with logo:







 **orient Honest**
MEDICAL INTERNATIONAL

