

### **The property of Disposable cotton soft towel**

**All cotton towel are made of the highest quality materials, providing excellent protection and superior comfort:**

- 1.Size: 20 x 20 cm each piece or customized
2. Made of cotton material, soft and comfortable, clean and sanitary.
3. Convenient and lightweight , easy to carry out .
4. Can be wet and dry, are great for fragile or sensitive skin
5. Great water absorption, skin-friendly and comfortable.
6. Cotton tissue can instead of face towel and will be much cleaner to dry your face after washing.

### **The Specification of Disposable cotton soft towel**

Name	Disposable cotton soft towel
Material	100% cotton;spunlace nonwoven
Style	Plain, mesh, pearl and customization patterns available
Size	15*20cm,20*20cm ,or other size can be customized
Color	white
Packing	1roll/portable bag ;100pcs/pack ,can be customized
Carton Box	All inner bag and carton printing can be as your design
Carton Size	52x38x30cm /52x19x30cm
Applications	personal care, makeup removing, house cleaning, travel use, camping, similar outdoor activities and business occasions.
MOQ	1000 bags
Sample	Can be supplied free for your quality checking within 3 working days
Payment way	T/T,PAYPAL,Westunion,L/C, D/A,D/P,ESCROW
Delivery time	as you quantity
FOB Port	Wuhan/Shanghai/Ningbo/Guangzhou

### **The picture of Disposable cotton soft towel**







orient Honest  
MEDICAL INTERNATIONAL



## Wet use

Wash Face  
Soft, Clean, No scent, No debris





orient Honest  
MEDICAL INTERNATIONAL



**Lock Abundant Liquid**