

The property of Disposable cotton tissue roll

All cotton towel are made of the highest quality materials, providing excellent protection and superior comfort:

- 1.Size: 20 x 20 cm each piece or customized
2. Made of cotton material, soft and comfortable, clean and sanitary.
3. Convenient and lightweight , easy to carry out .
4. Can be wet and dry, are great for fragile or sensitive skin
5. Great water absorption, skin-friendly and comfortable.
6. Cotton tissue can instead of face towel and will be much cleaner to dry your face after washing.

The Specification of Disposable cotton tissue roll

Name	Disposable cotton tissue roll
Material	100% cotton;spunlace nonwoven
Style	Plain, mesh, pearl and customization patterns available
Size	15*20cm,20*20cm ,or other size can be customized
Color	white
Packing	1roll/portable bag ;100pcs/pack ,can be customized
Carton Box	All inner bag and carton printing can be as your design
Carton Size	52x38x30cm /52x19x30cm
Applications	personal care, makeup removing, house cleaning, travel use, camping, similar outdoor activities and business occasions.
MOQ	1000 bags
Sample	Can be supplied free for your quality checking within 3 working days
Payment way	T/T,PAYPAL,Westunion,L/C, D/A,D/P,ESCROW
Delivery time	as you quantity
FOB Port	Wuhan/Shanghai/Ningbo/Guangzhou

The picture of Disposable cotton tissue roll



How to Use 1



Step 1

Remove the roll tube from the center of the towel.



Step 2

Pick up the first count of the towel from the center.



Step 3

Put it back in the original packing bag and help yourself.

How to Use 2



Step 1

Take the towel out of the original package.



Step 2

Get a holder or a rack, and place the product.



Step 3

Use it like you would use a kitchen towel.



