

Product details

The material of cute medical mask

- 1) 3layers: non woven spunbond + melt blown filter paper +spunbond nonwoven
- 2) The fabric has maximum fluid resistance and protection against body fluids
- 3) Bacterial Filtration Efficiency: 99.7%, Particle Filtration Efficiency: 98.5%,
Particle size: 0.1 μm

The performance of children face mask

- 1) Ultra-sonic weld ensures strongest constructed mask, eliminates holes and defects
- 2) Adjustable Nose Piece for better fit, very comfortable to wear

The trade terms of custom printed surgical mask

MOQ	500,000 pcs
Payment terms	T/T, L/C, D/A, D/P, Western Union, ESCROW,Paypal
Supply ability	1,000,000 pcs/day
Delivery time	According to the quantity
Packing details	50 pcs/box,2000 pcs/ctn, or as your requirement

Product display




orient Honest
MEDICAL INTERNATIONAL





orient Honest
MEDICAL INTERNATIONAL





The regular shipment way

