

## **3ply disposable surgical face mask , antibacterial face mask 17.5\*9.5cm , PP nonwoven soft anti-odor face mask**

### **Indications of 3ply disposable surgical face mask**

- 1) Sterile face mask is for hospital, clinical, medical personnel and patients;
- 2) Non-sterile face mask can be used for dust-free workshop, food factory and so on.

### **Specifications of antibacterial face mask**

Material	PP non woven + high filtration paper + PP non woven
Fabric weight	Could be as your request
Ear loop	Tie on or ear loop
Size	17.5*9.5cm, or as your requirements
Color	white/blue/green/pink/yellow/orange/purple /lemon/Capri blue/burgundy/black
Packing	50 pcs/box, 2000 pcs/ctn, or as your requirement
Packing design	all inner box and outer carton printing can be as your design
Applications	hospital, laboratory, food industry, dust-free workshop, electronic manufactures and so on

### **Packaging of anti-odor face mask**

50pcs/box, 40boxes/ctn, 52x19x42cm;  
1\*20 GP container can hold 1,400,000 Pcs of masks, 700 Ctns;  
1\*40 HC container can hold 3,450,000 Pcs of mask, 1725 Ctns.  
Or as your requirements.

### **Picture of 3ply disposable surgical face mask**

