

## Application of surgery mask

- 1.For hospital clinical medical and patients,dust free work shop,food factory.
- 2.To prevent of spread of body fluids and blood splash content.

## Parameters of surgical disposable face mask

- 1) Material: PP non woven + high filtration paper + PP non woven
- 2) Color: White/blue/green/pink/yellow/orange/purple /lemon/Capri blue/burgundy
- 3) Ear loop: tie on or elastic surgical disposable face mask
- 4) Size: 17.5\*9.5 cm for adult; 14.5\*9.5 for children; 12.5\*7.5cm for baby
- 5) Nose bar: plastic or metal, flexible,fitting fixed

## Packaging of medical face masks with design

Specification	Normal packing	Carton size	1*40 HQ container can hold (Pcs)
17.5*9.5cm,3ply, 18+20+20G	50 pcs/box,40 boxes/ctn	52*19*42cm	3,280,000 Pcs
17.5*9.5cm,3ply, 25+25+25G	50 pcs/box,40 boxes/ctn	51*37*30cm	2,400,000 Pcs
Remark: all can be special customized according to your actual requirements			

## Our factory

- 1) Manufacturer with 12 years history in producing surgical disposable face mask, and all kinds of disposable protective products;
- 2) Strong production power:1\*40 HQ container can be finished in 15 days;
- 3) Promising sales team: service in 24 hours
- 4) High quality insurance: ISO9001, ISO13485, CE certificated, FDA Registered.

## The trade terms of medical face masks with design

MOQ	100,000 pcs
Price	Negotiable
Payment terms	T/T, L/C, D/A, D/P, Western Union, ESCROW,Paypal
Supply ability	1,000,000 pcs/day,1*40 HQ container can be finished in 15 days

Delivery time	Normally 30 days after confirmed the order
Sample	Sample is free, can be prepared well in 2 days

**Picture of medical face masks with design**





orient Honest  
MEDICAL INTERNATIONAL

