

Application of surgery mask

1. For hospital clinical medical and patients, dust free work shop, food factory.
2. To prevent spread of body fluids and blood splash content.

Parameters of surgical disposable face mask

- 1) Material: PP non woven + high filtration paper + PP non woven
- 2) Color: White/blue/green/pink/yellow/orange/purple /lemon/Capri blue/burgundy
- 3) Ear loop: tie on or elastic surgical disposable face mask
- 4) Size: 17.5*9.5 cm for adult; 14.5*9.5 for children; 12.5*7.5cm for baby
- 5) Nose bar: plastic or metal, flexible, fitting fixed

Packaging of medical face masks with design

Specification	Normal packing	Carton size	1*40 HQ container can hold (Pcs)
17.5*9.5cm, 3ply, 18+20+20G	50 pcs/box, 40 boxes/ctn	52*19*42cm	3,280,000 Pcs
17.5*9.5cm, 3ply, 25+25+25G	50 pcs/box, 40 boxes/ctn	51*37*30cm	2,400,000 Pcs
Remark: all can be special customized according to your actual requirements			

Our factory

- 1) Manufacturer with 12 years history in producing surgical disposable face mask, and all kinds of disposable protective products;
- 2) Strong production power: 1*40 HQ container can be finished in 15 days;
- 3) Promising sales team: service in 24 hours
- 4) High quality insurance: ISO9001, ISO13485, CE certificated, FDA Registered.

The trade terms of medical face masks with design

MOQ	100,000 pcs
Price	Negotiable
Payment terms	T/T, L/C, D/A, D/P, Western Union, ESCROW, Paypal
Supply ability	1,000,000 pcs/day, 1*40 HQ container can be finished in 15 days

Delivery time	Normally 30 days after confirmed the order
Sample	Sample is free, can be prepared well in 2 days

Picture of medical face masks with design





orient Honest
MEDICAL INTERNATIONAL

