

## **Anti-bacteria surgical facial mask , cleanroom medical masks , non woven face mask supplier**

### **The feature of surgical facial mask**

- 1) Three-dimensional design: specified Anti-SARS masks,sterilization,anti-flu,dust
- 2) With nose clip: flexible,fitting fixed
- 3) Fine edge: the edge of full heat sealing technology package
- 4) Tie on ear loop: comfortable and fitted

### **Specifications of medical masks**

Material	PP non woven + high filtration paper + PP non woven
Fabric weight	Could be as your request
Ear loop	Tie on or ear loop medical masks
Size	17.5*9.5cm, or as your requirements
Color	white/blue/green/pink/yellow/orange/purple /lemon/Capri blue/burgundy/black
Packing	50 pcs/box,2000 pcs/ctn, or as your requirement
Packing design	all inner box and outer carton printing can be as your design
Applications	hospital,laboratory,food industry,dust-free workshop,electronic manufactures and so on

### **Payment terms of surgical facial mask**

MOQ	100,000 pcs
Price	Negotiable
Payment terms	T/T, L/C, D/A, D/P, Western Union, ESCROW,Paypal
Supply ability	1,000,000 pcs/day,1*40 HQ container can be finished in 15 days
Delivery time	Normally 30 days after confirmed the order
Sample	Sample is free, can be prepared well in 2 days

### **Picture of non woven face mask**

