

Anti-bacteria surgical type mask , polypropylene non woven mask , elastic cord hospital mask manufacturer

Overview of surgical type mask

- 1) Looped Ears For Secure Hold
- 2) One size to fit all, approximate size: 17.5 cm x 9.5 cm
- 3) Covers Nose, Mouth And Chin For Full Protection
- 4) Color: white/blue/green/pink/yellow/orange/purple /lemon/Capri blue/burgundy/black, we have 10 colors totally
- 5) NELSON,ISO,CE certificates, BFE>99%,type II /Type IIR standard
- 6) Latex free,easy breath,eco friendly
- 7) Multi usage: surgical, flu, dust protectors, dentists, salons, general maintenance etc

Packaging of non woven mask

50pcs/box,40boxes/ctn,52x38x30cm;
1*20 GP container can hold 1,000,000 Pcs of masks;
1*40 HC container can hold 2,450,000 Pcs of mask.

We can also do packaging as your requirements.

The trade terms of hospital mask

MOQ	100,000 pcs
Price	Negotiable
Payment terms	T/T, L/C, D/A, D/P, Western Union, ESCROW,Paypal
Supply ability	1,000,000 pcs/day,1*40 HQ container can be finished in 15 days
Delivery time	Normally 30 days after confirmed the order
Sample	Sample is free, can be prepared well in 2 days

Picture of hospital mask



orient Honest
MEDICAL INTERNATIONAL

