

Indications of CPE anti skid shoe cover:

CPE medical disposable shoe cover for hospital , 45x16cm , blue can be used for hospital, dust-free workshop, laboratory, food industry, electronic manufactures and so on.

It also could be used for hospital clinical medical personnel and patients.

Product Advantages:

- 1, Available material: CPE
- 2, Size: 45x16cm ,we also can make 38x15cm, 40x15cm, 40x17cm, 45*18cm
- 3, Good protective, lightweight, latex free
- 4, Use sewing or automatic technology, comfortable wearing
- 5, Available color: blue ,white, green, yellow, red

Quality Control:

1. We will not start to make the products until you confirm the sample.
2. The packing will be water-proof, damp-proof and sealed

Specification of anti slip CPE shoe cover:

Item Name	CPE shoe cover
Type:	Machine-made or sewing
Material	CPE1.5-6g
Color	Blue, white, etc.
Size:	38x15cm,40x15cm,41x17cm.....
Packing:	100pcs/bag, 2000pcs/ctn
Applications	Household, hotel, electronics factory, pharmaceutical factory , hospital and other cleaning work environment
Daily output	80,000 pcs/day
Payment terms	T/T, PayPal, Western union, L/C,D/A,D/P,ESCROW,
Delivery time	Capacity: 200,000pcs/day, 1*40HQ can be finished within 10-15 days;
Certificate	CE / ISO13485
FOB port	Wuhan/Shanghai/Ningbo/Guangzhou

Details picture of CPE shoe cover:



