

All shoe cover are made of the highest quality materials, providing excellent protection and superior comfort:

The Specification of [shoe cover](#)

Increase the size of shoe cover by 15*40 to meet the size demand of export foot
The thickening of shoe covers with 3.0G new material meets the needs of customers

Product name: disposable [CPE plastic shoe cover](#)

Grams of products: 2.5g,3.0g,3.5g,4.0g and other multi-specification orders

Product packaging: 100 PCS/bag 4000 PCS/box If you want to calculate the size of the space box, please contact customer service

Product style: CPE plastic /PE plastic

Product Description	
Item name	Disposable PE shoe covers
Size	15X36CM,15X38CM,15X40CM,17X45CM,Can be Customed
Material	PE
Application	Hospital, home, computer room, clean room, laboratory
Remark	Auto Machine Made/Hand Made
Colour	Blue,red or as request
Package	100pca/bag,4000pcs/ctn
Sample	Free of charge
MOQ	100000PCS
Remark	All can be special customized according to your actual





FAQ

1.How can I get some samples?

Low sample charge

2.What is the sampling time?

Usually we need 7 days to prepare the samples.Sometimes we make samples according to clients' requirements.

3.Can I print our LOGO on package?

Sure,customized package and OEM printing are available.

4.How about the payment for mass production?

Normally we take 50% earnest in advance,and take final payment before delivery.

5.How do you control the quality of your products?

We always do our best to satisfy our clients' requirements. Because we know the most important thing is to meet the market's need.So we will do our best all the time. From raw material to final packaging, our quality control department will responsible for quality checking in each process.