

China best supplier for hospital surgical beard mask , disposable beard face masks , poly beard cover

Indications of surgical beard mask

It's widely used in hospital, hotel, restaurant, food process, catering and industrial sectors.

Fabric material of disposable beard face masks

100% polypropylene non-woven, non-irritating and soft

Size: 19", 21"

Packaging of beard cover

100 Pcs/PE bag , 1000 Pcs/carton, Ctn size: 26*22*32cm

100 Pcs/PE bag , 2000 Pcs/carton, Ctn size: 52*23*26cm

Storage of surgical beard mask

- 1) The product will maintain its properties when it is stored in a dry place at temperatures between 10 ° and 30 ° C.
- 2) Not exposed to direct sunlight.

Discarded

- 1) Once finished use can be recycled or incinerated.
- 2) Packaging should be deposited in the appropriate container for recycling

Trade terms of disposable beard face masks

MOQ	100,000 pcs
Price	Negotiable
Payment terms	T/T, L/C, D/A, D/P, Western Union, ESCROW, Paypal
Supply ability	1,000,000 pcs/day, 1*40 HQ container can be finished in 15 days
Delivery time	Normally 30 days after confirmed the order
Sample	Sample is free, can be prepared well in 2 days

Picture of beard cover

