

Introduction of disposable medical dental [cotton roll](#):

1. high absorbing capacity
2. made from 100% high-purity absorbent cotton, pure white
3. Packaged in plastic for protection
4. Size:
 - 1# 8mm X 38mm (0.315"X1-1/2")
 - 2# 10mm X 38mm (0.375"X1-1/2")
 - 3# 10mm X 38mm (0.375"X1-1/2")
 - 4# 10mm X 38mm (0.375"X1-1/2")
 - 5# 12mm X 38mm (0.472"X1-1/2")
 - 6# 15mm X 38mm (0.551"X1-1/2")

Advantages of [cotton roll](#)

1. Soft, pliable, no cellulose rayon fibers, non-linting, and pleasant for patient
2. absorbing blood and exudates or for cleaning wound.
3. High absorbency and firmly attached
4. Will not adhere to mucous membranes
5. Maintain shape better when wet

Application of [cotton roll](#):

For dental clinic or hospital.

Packing:

19 pcs/roll, 36 roll/bag or customized

Specification of dental cotton roll:

Material	100% high-purity absorbent cotton
Feature	dental cotton roll
Color	pure white
Size	1#: 8*38mm, 2#: 10*38mm, 3#: 12*38mm, 4#: 15*38mm
Delivery	within 15 days after receipt of deposit.

Terms of payment	T/T, L/C, Western Union, Escrow, Paypal, etc.
Storage	Stored in dry, humidity below 80%, ventilated, non-corrosive gases warehouse
Note	1.This product is Non-sterile 2.Validity: 2 years 3.Disposable
Certification	CE0197, ISO9001, ISO13485
OEM	1.Material or other specifications can be according to customers' requirements. 2.Customized Logo/brand printed. 3.Customized packaging available.





orient Honest
MEDICAL INTERNATIONAL





orient Honest
MEDICAL INTERNATIONAL

