

Introduction of disposable medical dental [cotton roll](#):

- 1.high absorbing capacity
- 2.made from 100% high-purity absorbent cotton, pure white
- 3.Packaged in plastic for protection
- 4.Size:
 - 1 1# 8mm X 38mm (0.315"X1-1/2")
 - 2 2# 10mm X 38mm(0.375"X1-1/2")
 - 3 3# 10mm X38mm(0.375"X1-1/2")
 - 4 4# 10mm X38mm(0.375"X1-1/2")
 - 5 5# 12mm X 38mm(0.472"X1-1/2")
 - 6 6# 15mm X 38mm(0.551"X1-1/2")

Advantages of [cotton roll](#)

- 1.Soft,pliable,no cellulose rayon fibers, non-linting, and pleasant for patient
- 2.absorbing blood and exudates or for cleaning wound.
- 3.High absorbency and firmly attached
- 4.Will not adhere to mucous membranes
5. Maintain shape better when wet

Application of [cotton roll](#):

For dental clinic or hospital.

Packing:

19pcs/roll,36roll/bag or customized

Specification of dental cotton roll:

Material	100% high-purity absorbent cotton
Feature	dental cotton roll
Color	pure white
Size	1#:8*38mm,2#:10*38mm,3#:12*38mm,4#:15*38mm
Delivery	within 15 days after receipt of deposit.

Terms of payment	T/T, L/C, Western Union, Escrow, Paypal, etc.
Storage	Stored in dry, humidity below 80%, ventilated, non-corrosive gases warehouse
Note	1.This product is Non-sterile 2.Validity: 2 years 3.Disposable
Certification	CE0197, ISO9001, ISO13485
OEM	1.Material or other specifications can be according to customers' requirements. 2.Customized Logo/brand printed. 3.Customized packaging available.





orient Honest
MEDICAL INTERNATIONAL





orient Honest
MEDICAL INTERNATIONAL

