

China disposable non woven surgical breathing ear loop medical mask suppliers

Normal packing

- 1.Non-sterile: 50pcs/box,40boxes/ctn,52x38x30cm.
1*20FT can load:500 cases(1000,000 pcs),
1*40HQ can load:1150 cases(2300,000 pcs).
- 2.Sterile: 1pcs/bag, 50bags/box,1000pcs/ctn, carton size:52x19x30cm.

Indication

- 1.Sterile products are used for hospital clinical medical personnel and patients.
- 2.Non-sterile products can be used for dust-free workshop,hospital, laboratory,
- 3.To prevent the spread of body fluids and blood splash content and isolate dust,

Specifications

Product Name	China disposable non woven surgical breathing ear loop medical mask suppliers
Model No	OH-2001
Material	3ply, PP non woven + Filter paper + PP non woven
Size	17.5x9.5cm(adult)
Style	Elastic ear-loop(latex free)
Layer	3ply
Color	white
Packing	50pcs/box,40boxes/ctn,2000pcs/ctn. 1*20FT can load:500 cases(1000,000 pcs), 1*40HQ can load:1150 cases(2300,000 pcs).
Carton Size	52*38*30cm
GW	GW: 7.5kg/ctn, NW: 7kg/ctn.
MOQ	5000PCS
Certificate	CE / ISO13485 / NELSON / FDA

Photos





orient Honest
MEDICAL INTERNATIONAL



