

## **Introduction of disposable [cotton ball](#):**

1. Gauze balls and Non-woven ball are ideal selections for the absorption of blood and exudate.
2. with or without elastic ring.
3. We can provide the one with X ray or without X ray.
4. It can also be used for wound cleansing and for skin disinfection.

## **[Cotton ball](#):**

It is made of cotton

## **Advantages**

1. 100% cotton, pure white, absorbent & soft
2. CE and ISO certified
3. Factory direct

## **Indications:**

For hospital clinical medical and patients.

## **Packing:**

100pcs/bag, 2000pcs/carton, or as customer request

## **Specification of [cotton ball](#):**

Properties:	Medical Materials & Accessories
Type:	Surgical Supplies
Model Number:	Request
Material:	100% Cotton
Product name:	Gauze Balls
Color:	White
Cotton yarn:	21's, 32's, 40's
mesh:	26x18, 30x20, 30x30, Etc.
Packing:	100, 250, 500/polybag, 3, 5, 10, 20, 30pcs/pouch
Certificate:	CE/ISO13485
Usage:	For wound cleansing and for skin disinfection





orient Honest  
MEDICAL INTERNATIONAL





orient Honest  
MEDICAL INTERNATIONAL





orient Honest  
MEDICAL INTERNATIONAL







orient Honest  
MEDICAL INTERNATIONAL

