

Product details

Introduction of China exporter for green SMS lab coats:

1. Material: SMS
2. Color: green, etc.
3. Weight: 45gsm
4. Style:knitted Collar and cuffs, with pockets
5. Size:S-2XL.

Specifications of China exporter for green SMS lab coats:

Name	China exporter for green SMS lab coats,three pockets disposable hospital lab coat,hot selling SMS lab coats wholesale
Material grammage	SMS45g
Style	with knit collar and cuffs, with three pockets
Size	S M L XL 2XL
Color	green, or blue, yellow, pink, red, grey, black etc.
Packing	Non sterile:5pcs/bag,10 bags/ctnor ; Sterile:1pcs/bag, 50pcs/ctn
Packing design	All inner box and carton printing can be as your design
MOQ	5,000 pcs
Sample	Can be supplied free for your quality checking within 3 days
Payment way	T/T, PayPal, Western union, L/C,D/A,D/P,ESCROW,
Delivery time	Capacity: 10,000pcs/day, 1*40HQ can be finished within 10-15 days;
FOB port	Wuhan/Shanghai/Ningbo/Guangzhou

Features of China exporter for green SMS lab coats:

SMS visitor coat is made from hydrophobic SMS material ,Latex-free; abrasion-resistant; low lint; with a high level of fluid repellency; a good barrier to blood, body fluids and pathogens.

Product display



orient Honest
MEDICAL INTERNATIONAL






orient Honest
MEDICAL INTERNATIONAL



Our business process



The regular shipment way



FAQ

Q: 1. Are you a manufacturer or a trading company, and when was your company established?

A: Yes, we are directly manufacturer, our factory was established in 2003, we have our own exported rights company in 2013.

Q: 2. What is the annual turnover for your company?

A: Normally it is around USD18 million-20 million, it is still in the growing.

Q: 3. Do your company have any quality management system such as ISO, EN or CE certificate(s)?

A: Yes, we have CE, ISO13485, FDA, REACH and NELSON report.

Q: 4. Can I get the free samples for quality checking? How long time can I get it?

A: Yes, our samples can be supplied free within 3-7 days.