

Product details

Features of China factory for lab coat with knit cuffs:

1. Material: PP50g
2. Color: green
3. Grammage: 50gsm
4. Style: with knitted cuffs/collar, with two pockets
5. Size:S-2XL.
6. Lightweight, durable and superior breathable; Allow great freedom of movement.

Specifications of China factory for lab coat with knit cuffs:

Name	China factory for lab coat with knit cuffs□disposable green lab coat with pockets□nonwoven PP50g lab coats wholesale for adults
Material	PP50g
Style	With two pockets, knitted collar, button on the front, knitted cuff or other styles
Size	S M L XL 2XL,3XL,etc
Color	green, etc.
Packing	1-5pcs/bag,10 bags per carton
Packing design	All inner box and carton printing can be as your design
MOQ	5,000 pcs
Sample	Can be supplied free for your quality checking within 3 days
Payment way	T/T, PayPal, Western union, L/C,D/A,D/P,ESCROW,
Delivery time	Capacity: 10,000pcs/day, 1*40HQ can be finished within 10-15 days;
FOB port	Wuhan/Shanghai/Ningbo/Guangzhou

Product display



orient Honest
MEDICAL INTERNATIONAL






orient Honest
MEDICAL INTERNATIONAL



Our business process



The regular shipment way



FAQ

Q: 1. Are you a manufacturer or a trading company, and when was your company established?

A: Yes, we are directly manufacturer, our factory was established in 2003, we have our own exported rights company in 2013.

Q: 2. What is the annual turnover for your company?

A: Normally it is around USD18 million-20 million, it is still in the growing.

Q: 3. Do your company have any quality management system such as ISO, EN or CE certificate(s)?

A: Yes, we have CE, ISO13485, FDA. REACH and NELSON report.

Q: 4.Can I get the free samples for quality checking? How long time can I get it?

A: Yes, our samples can be supplied free within 3-7 days.