

Product details

China manufacturer disposable operating gowns:

China manufacturer disposable operating gowns ,disposable non-sterile surgical gown,PP+PE disposable surgical gown is made from PP+PE material and designed for optimal protection during surgical procedures of moderate duration involving moderate amounts of body fluids and moderate risk of infection.

Specification of China manufacturer disposable operating gowns:

material	PP+PE, other disposable materials also could be offered
Weight	PP+PE 45 gsm
Style	tie/velcro on collar,tie on waist,elastic/knitted cuff
Size	120x130cm,130x150cm,150x160cm,150x190cm,...or could be customized
Color	green, etc.
packing	non sterile,10 pcs/bag,10 bags/ctn. sterile,1 pcs/sterile bag,100 pcs/ctn
packing design	all inner box and carton printing can be as your design
carton size	58X28X48 cm
Applications	hospital clinical medical personnel and patients. this kind of surgical gowns also could be used in dust-free workshop,laboratory, food industry, electronic manufactures and so on
MOQ	5000 pcs
Sample	can be supplied free for your quality checking within 3 days

Product display



orient Honest
MEDICAL INTERNATIONAL






orient Honest
MEDICAL INTERNATIONAL



Our business process



The regular shipment way



FAQ

Q: 1. Are you a manufacturer or a trading company, and when was your company established?

A: Yes, we are directly manufacturer, our factory was established in 2003, we have our own exported rights company in 2013.

Q: 2. What is the annual turnover for your company?

A: Normally it is around USD18 million-20 million, it is still in the growing.

Q: 3. Do your company have any quality management system such as ISO, EN or CE certificate(s)?

A: Yes, we have CE, ISO13485, FDA, REACH and NELSON report.

Q: 4. Can I get the free samples for quality checking? How long time can I get it?

A: Yes, our samples can be supplied free within 3-7 days.

