

## Product details

### Introduction of China manufacturer non-woven isolation gown

Disposable surgical gown is made from hydrophobic polypropylene material with PE coated, or SMS/Spunlace material with PE coated.

### The packaging of China manufacturer non-woven isolation gown:

Sterile: 1 pcs/paper plastic bag,50 pcs/Ctn

Non-sterilize: 5 pcs/bag, 50 pcs/ctn

Or as your requirements.

### Specifications of China manufacturer non-woven isolation gown

Name:	China manufacturer non-woven isolation gown,large size doctor nonwoven surgical gown, one time use non-woven surgical gown
Material:	SMS
Grammage:	45 gsm
Closure:	with belt closure, or with Velcro closing
Application:	for medical, surgical, pharmaceutical, laboratorial, industrial sector
Size:	S/M/L/XL/XXL/XXXL,
Color:	Green,Blue, etc, optional
Packing:	Normal packing:1pcs/polybag,50pcs/carton
Storage:	Stored in dry, humidity below 80%, ventilated, non-corrosive gases warehouse
Certification:	CE, ISO13485,FDA

## Product display




orient Honest  
MEDICAL INTERNATIONAL






orient Honest  
MEDICAL INTERNATIONAL





## Our business process



## The regular shipment way



## FAQ

- 1, Free sample.
- 2, Small order acceptable;
- 3, L/C, T/T, D/P terms available;
4. FOB, CNF, CIF, EXW terms available;
- 5, Competitive and stable price because we are factory.
- 6, Qualified products. 3rd part product inspection is acceptable.
- 7, OEM service. More options on material, pattern, size and packing;
- 8, Timely delivery. The best lead time is 25-30 day for a 40HQ container.
- 9, Professional export service. Documentation, clearance, sea transportation, all in one.
- 10, Good after sale service. We take care complaints as we would like to be your long-term supplier!