

Introduction of China manufacturer surgical doctors caps:

Certification: CE / ISO13485 / FDA

Place of Origin: China

Material : PP nonwoven

Color: green,blue,pink,yellow,pink,orange,purple,red, white,(all colors can be available)

Style: double elastic clip caps

Weight: 10g, or 11g-30g

Size: 21''

Packing details: 100 pcs/bag,20 pcs/ctn.

Features of China manufacturer surgical doctors caps:

1. Ecofriendly, breathable, health, nonirritant, economical, Soft, breathable, lightweight and comfor table .
2. Different design of products is available on your request.
3. Samples and customization available.
4. Latex free, class one flammability rating.

Specifications of China manufacturer surgical doctors caps:

| | |
|----------|---|
| Material | PP non woven |
| Weight | 10-30g |
| Style | Elastic (latex free), single and double elastic style |
| Size | 18'',19'',21'',24''28'' |

| | |
|----------------|--|
| Color | blue,pink,yellow,pink,orange,purple,red,white,green,(all colors can be available) |
| packing | 100pcs/bag, 20bags/ctn, 2000pcs/ctn |
| packing design | All inner box and carton printing can be as your design |
| carton size | 58x25x27cm |
| Applications | hospital clinical medical personnel and patients. dust-free workshop,, laboratory, food industry, electronic manufactures and so on |
| MOQ | 100,000 PCS |

Pictures of China manufacturer surgical doctors caps:



orient Honest
MEDICAL INTERNATIONAL



