

**Introduction of China manufacturer surgical doctors caps:**

Certification: CE / ISO13485 / FDA

Place of Origin: China

Material : PP nonwoven

Color: green,blue,pink,yellow,pink,orange,purple,red, white,(all colors can be available)

Style: double elastic clip caps

Weight: 10g, or 11g-30g

Size: 21''

Packing details: 100 pcs/bag,20 pcs/ctn.

**Features of China manufacturer surgical doctors caps:**

- 1. Ecofriendly, breathable, health, nonirritant, economical, Soft, breathable, lightweight and comfor table .
- 2. Different design of products is available on your request.
- 3. Samples and customization available.
- 4. Latex free, class one flammability rating.

**Specifications of China manufacturer surgical doctors caps:**

Material	PP non woven
Weight	10-30g
Style	Elastic (latex free), single and double elastic style
Size	18",19",21",24"28"

Color	blue,pink,yellow,pink,orange,purple,red,white,green,(all colors can be available)
packing	100pcs/bag, 20bags/ctn, 2000pcs/ctn
packing design	All inner box and carton printing can be as your design
carton size	58x25x27cm
Applications	hospital clinical medical personnel and patients. dust-free workshop,, laboratory, food industry, electronic manufactures and so on
MOQ	100,000 PCS

### **Pictures of China manufacturer surgical doctors caps:**



