

## Introduction of China plant for 19" single elastic clip cap:

China plant for 19" single elastic clip cap made by SBPP nonwoven material and pleated elastic band , certified CE and ISO13485. Also called disposable strip cap. Use for the dentist and doctor to protect the hair clean for their daily working.

## Specifications of China plant for 19" single elastic clip cap:

<b>Product:</b>	China plant for 19" single elastic clip cap , nonwoven dental clip caps, dentist clip mob bouffant cap
<b>Material:</b>	Nonwoven pp with elastic
<b>Grammage:</b>	10g~16g
<b>Style:</b>	Single elastic or double elastic
<b>Size:</b>	19", or 18", 20", 21", etc.
<b>Color:</b>	White, Green, yellow, Blue, etc.
<b>MOQ</b>	100,000pcs
<b>Packaging:</b>	100pcs/bag ; 2000pcs/ctn
<b>Application:</b>	for household, medical, surgical, pharmaceutical, laboratorial, industrial sectors, painting, cleaning room, etc.
<b>Storage:</b>	Stored in dry, humidity below 80%, ventilated, non-corrosive gases warehouse
<b>Terms of payment:</b>	T/T, L/C, Paypal, D/P, D/A, etc.
<b>Other service:</b>	OEM and Customized Logo/brand
<b>Certification:</b>	CE, ISO13485

## Application:

China plant for 19" single elastic clip cap , nonwoven dental clip caps, dentist clip mob bouffant cap can be used for household, hospital, laboratory, foods, electronic industry. Prevent and isolate dust, particles, etc., can be a visit protection. Economical and disposable. Conform to the hygienic standard.

### **Packing:**

100 pcs/bag,20 bags/ctn,carton size:58x25x27cm.

1\*20 container can load 715ctns,total 1430,000 pcs.

1\*40 container can load 1737ctns,total 3474,000 pcs.

### **Sample Policy:**

1. Sample according to your design drawing. Sampling time:3-5days .
2. Design your own sample, then remake it.
3. Existing Samples Sampling time:1day.
4. Make mass production after you confirm the sample.

### **Pictures of China plant for 19" single elastic clip cap:**





orient Honest  
MEDICAL INTERNATIONAL



