

Introduction of China plant for 19" single elastic clip cap:

China plant for 19" single elastic clip cap made by SBPP nonwoven material and pleated elastic band , certified CE and ISO13485. Also called disposable strip cap. Use for the dentist and doctor to protect the hair clean for their daily working.

Specifications of China plant for 19" single elastic clip cap:

Product:

China plant for 19" single elastic clip cap , nonwoven dental clip caps, dentist clip mob bouffant cap

Material:

Nonwoven pp with elastic

Grammage:

10g~16g

Style:

Single elastic or double elastic

Size:

19", or 18", 20", 21", etc.

Color:

White, Green, yellow, Blue, etc.

MOQ

100,000pcs

Packaging:

100pcs/bag ; 2000pcs/ctn

Application:

for household, medical, surgical, pharmaceutical, laboratorial, industrial sectors, painting, cleaning room, etc.

Storage:

Stored in dry, humidity below 80%, ventilated, non-corrosive gases warehouse

Terms of payment:

T/T, L/C, Paypal, D/P, D/A, etc.

Other service:

OEM and Customized Logo/brand

Certification:

CE, ISO13485

Application:

China plant for 19" single elastic clip cap , nonwoven dental clip caps, dentist clip mob bouffant cap can be used for household, hospital, laboratory, foods, electronic industry. Prevent and isolate dust, particles, etc., can be a visit protection. Economical and disposable. Conform to the hygienic standard.

Packing:

100 pcs/bag,20 bags/ctn,carton size:58x25x27cm.

1*20 container can load 715ctns,total 1430,000 pcs.

1*40 container can load 1737ctns,total 3474,000 pcs.

Sample Policy:

1. Sample according to your design drawing. Sampling time:3-5days .
2. Design your own sample, then remake it.
3. Existing Samples Sampling time:1day.
4. Make mass production after you confirm the sample.

Pictures of China plant for 19" single elastic clip cap:





