

Product details

Description of China workshop food industry lab coat

1. Material:pp40g
2. Color: White,
3. Style: knitted cuff,knit collar,with pockets
4. Size: S-XXXL sizes available
5. Packing: 5PCS/polybag, 50PCS/carton and as you require
- 6.1*20" FCL: 50000PCS (tolerance as different sizes and materials)
7. Mini order quantity: 5000PCS
8. Certificate: CE, FDA, ISO13485

Specifications of China workshop food industry lab coat

Material	pp non woven
Grammage	40g ,Could as your request
Style	With knit collar, with knitted cuffs,with pockets
Size	S-6XL
Color	white□ etc.
Packing	5pcs/box,50 bags/ctn, or as your requirement
Packing design	all inner bag and carton printing can be as your design
Applications	dustfree workshop, hospital, laboratory, food industry, electronic manufactures and so on

Product display



orient Honest
MEDICAL INTERNATIONAL






orient Honest
MEDICAL INTERNATIONAL



Our business process



The regular shipment way



FAQ

Q: 1.Are you a manufacturer or a trading company, and when was your company established?

A: Yes, we are directly manufacturer, our factory was established in 2003, we have our own exported rights company in 2013.

Q: 2.What is the annual turnover for your company?

A: Normally it is around USD18 million-20 million, it is still in the growing.

Q: 3. Do your company have any quality management system such as ISO, EN or CE certificate(s)?

A: Yes, we have CE,ISO13485,FDA. REACH and NELSON report.

Q: 4.Can I get the free samples for quality checking? How long time can I get it?

A: Yes, our samples can be supplied free within 3-7 days.